




Patientinstruktion för rektal provtagning klamydia och gonorré -PCR



1. För in bomullspinnen 2-3 cm i analöppningen, precis förbi slutmuskeln. Roter pinnen mot slemhinnans vägg i 10-30 sekunder. Om det finns synlig avföring på pinnen ska du ta nytt prov med en ny provtagningspinne.
2. Sätt pinnen i provtagningsröret med bomullsändan nedåt – bryt av skaftet vid den markerade skåran på pinnen.
3. Sätt på korken på röret och skruva åt ordentligt. Placera röret på anvisad plats



Patient instruction for rectal swab specimen collection for chlamydia and gonorrhoea -PCR



1. Insert the swab about 2-3 cm into the anal canal, just passed the anal sphincter muscle. Gently turn the swab for about 10-30 seconds while running the swab against the walls of the rectum. If the swab is grossly contaminated with faeces, discard and repeat the collection. Withdraw the swab carefully. Do not let the swab touch any surface before placing it into the collection tube.
2. Place the swab in the test tube with the cotton bud downwards – break the shaft of the swab at the marked point.
3. Screw on the lid of tube and make sure it is secure. Place the tube in the specified place.

